



COMPATIBLE WITH ALL OUR WATER & AIR SENSORS

INOVIEW V2


3 product ranges, 1 platform

Inoview is the **data visualization platform** for consulting your measurements in real time. Thanks to ongoing development, it now supports both our water and air monitoring devices – including the **Neptis** and **Aquaread** ranges for water and **Ostaro** for air quality. Hosted on a new, more secure and faster server, Inoview brings power and reliability to your environmental monitoring networks.

-  **REAL-TIME DATA VISUALIZATION**
-  **USER-FRIENDLY INTERFACE**
-  **ALARM AND ALERT MONITORING**
-  **EXPORT CUSTOMIZED REPORTS**
-  **CUSTOMIZABLE DASHBOARD**
-  **SECURE ACCESS AND DATA STORAGE**



HOME PAGE



GROUP INSTRUMENTS BY SITE
FOR EASIER PROJECT
MANAGEMENT


After logging in to your account using a secure login, you can view the data and status of your instruments in **real time** to ensure that your monitoring network is running smoothly:

- Map with **geolocation** of instruments
- Recent alerts sent by email (SMS optional) and triggered when user-defined parameter **thresholds** are exceeded
- List of the latest maintenance alarms issued on the **status** of instruments
- Latest **field notes** entered by operators

INSTRUMENT PAGE

Take a closer look at each of your instruments:

- Probe and modem **battery** status
- **GPS** coordinates and geolocation on map
- Last **value** received for each measured parameter
- Historic **graphs** and **tables** to zoom in on specific sections of data and see the value for each individual data point



COMPARISON PAGE

Compare your instruments or the values of various measured parameters for a fast, in-platform data analysis:

- Compare the same parameter across up to **6 instruments**
- Compare up to **6 parameters** from a single instrument



Let's talk!



SALES DEPARTMENT



+33 2 97 36 11 81



info.instrumentation@nke.fr